

DSLST protocol for *Pseudomonas aeruginosa*

A detailed protocol can be requested at contact@dlst.org.

The DLST is based on the analysis of 400bp and 350bp of the two *Pseudomonas aeruginosa* highly variable loci, *ms172* and *ms217*.

Chromosomal DNA extraction of *P. aeruginosa*

The DNA extraction protocol is based on a simple heating process in a 96-wells plate: bacterial isolates are incubated overnight at 37°C in 5ml of LB medium and crude DNA extracts are prepared by mixing 10µl of this broth with 90µl of TE and heating this solution at 95°C for 10 minutes. This reduces time and cost per DNA extraction. However, traditional DNA extraction methods for *P. aeruginosa* are also fully appropriate for DLST.

PCR amplification of the variable region of the *ms172* and *ms217* loci

This PCR protocol is used to amplify highly polymorphic regions of the *ms172* and *ms217* loci. The priming regions are localized on conserved regions outside of the polymorphic regions. The forward and reverse primers are:

<i>ms172</i> .forward (5'-3')	GGATTCTCTCGCACGAGGT
<i>ms172</i> .reverse (5'-3')	TACGTGACCTGACGTTGGTG
<i>ms217</i> .forward (5'-3')	TTCTGGCTGTCGCGACTGAT
<i>ms217</i> .reverse (5'-3')	GAACAGCGTCTTTTCCTCGC

PCR amplification of each loci (*ms172* and *ms217* separately) were carried out in 22µl reaction volumes containing 2µl of the DNA extract solution, 1U of *Taq* DNA polymerase (Invitrogen), 1X *Taq* Reaction Buffer, 1.5mM MgCl₂, 0.3µM of each primer, and 0.2mM of each dNTP. Cycling conditions were 5 min of initial denaturation at 94°C; 35 cycles consisting of 30s at 94°C, 30s at 60°C (*ms172*) or 64°C (*ms217*), 45s at 72°C; and a final extension for 10min at 72°C. Some optimization may be required according to the material used.

Sequencing of the reverse fragment

DNA fragments of 400bp for *ms172* and 350bp for *ms217* have been selected for the *P. aeruginosa* DLST typing method and the protocol has been adapted for the 3' sequencing (i.e. using the reverse primer) of these fragments. Annealing temperature for the sequencing reaction is 60 °C for the *ms172* and 64°C for the *ms217*.

Trimming patterns

Both sequences must be trimmed after a start signature that can be recognized on the beginning of the sequence:

Locus	start signature	length
ms172	ACGATGCTGGADCCA	400bp
ms217	CAGCATGGCG	350bp

For shorter alleles (sequence repeatedly shorter than 400bp or 350bp), please [contact](#) the curator of the website.

Allele and DLST type assignment

Once you obtained a sequence for both *ms172* and *ms217*, it is possible to obtain the allelic number corresponding to your sequence using the following options:

Submission of trace files.

A trace file can be submitted for the *ms172* and *ms217* loci. For both loci, the software will look for start signatures and will trim the sequences to their expected length (400 and 350bp). The sequences are then compared to the *ms172* and *ms217* database of alleles. If alleles that are already known are found, the corresponding alleles are given. . If one or two of the sequences are not known, the quality of the trace file is checked. A new allele number is given only if the quality of the trace file is good enough and if you have requested the permission to add new alleles using the [contact](#) form. If the quality is not good enough, a message indicate that a good quality trace file should be entered.

Submission of fasta files.

Fasta files with record named *ms172* and *ms217* for one isolate at a time can be submitted. The fasta file must be in the following form:

```
>ms172  
Sequence of 400bp  
>ms217  
Sequence of 350bp
```

It will soon be possible to submit fasta files with batches of isolates.

Submission of plain text sequences.

Plain text sequences can be entered for the *ms172* and *ms217* loci. These sequences must start with the corresponding start signature and must be trimmed at the correct length.